Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) Summary Data

Beijing Hotgen Biotech Co., Ltd.

Explanation of The Export License

To whom it may concern

According to the No.12 Government Notice published by MOC, GAC and NMPA, CCCMHPIE announced the list of companies which were permitted to export novel coronavirus diagnostic products, Hotgen is in the list as aforesaid.

The screenshot shows the website of the China Council for the Promotion of International Trade (CCPIT) China Pharmaceutical and Health Products Import and Export Association (CCCMHPIE). The page displays a table with one entry:

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

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Product Features

- High Accuracy, Specificity and Sensitivity ●
- No need instrument, get results in 15 minutes ●
- Room temperature storage ●
- Sample : Human Anterior Nares Swab ●
- Detect the presence of viral proteins ●
- Identify acute or early infection ●

Clinical Performance

(Disease Course 5-7 Days)

Sensitivity: 95.37%; Specificity: 99.13%; Accuracy: 97.31%.

Novel Coronavirus 2019-nCoV Antigen Test(Colloidal Gold)


Specimen Requirements

1 Sample collection

Gently insert the entire soft tip of the swab into one nostril for 1.5cm until you feel a bit of resistance.


Using medium pressure, rub the swab slowly in a circular motion around the inside wall of your nostril 4 times for a total time of 15 seconds.

Repeat the same process with the same swab in the other nostril.

2 Sample treatment

The swab after sampling is soaked below the liquid level of the sample extraction buffer. Rotate and press 3 times. The swab soaking time is not less than 15s.

The swab head is pressed, then take out the swab and tighten the sampling tube.

3 Sample preservation

The treated sample should be tested within 1h.

Test Procedure

Place the test cassette, sample extraction buffer at room temperature for 15-30 minutes, and equilibrate to room temperature (10-30 C).

Open the aluminum foil pouch of the test cassette, place the test cassette on a flat surface.


Add 4 drops of the treated sample into the sample well of the test cassette. (In case of chromatographic abnormalities, extra add 1-2 drops of the treated sample accordingly.) Incubate at 10-30 C for 15 minutes.


Observe the results after incubate at 10-30 C for 15 minutes. The result after 30 minutes is invalid.


Interpretation of result

Positive

Negative

Invalid

Clinical Performance

This study enrolled a total of 223 human anterior nasal swab specimens, of which 108 were positives for RT-PCR Test, 115 negatives for RT-PCR Test; As for data collection of the corresponding RT-PCR Test results.

Assessment system Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)	Reference system (clinical diagnostic results)		
	Positive(+)	Negative(-)	Total
Positive(+)	103	1	105
Negative(-)	5	114	118
Total	108	115	223

Sensitivity: 95.37%; Specificity: 99.13%; Accuracy: 97.31%.

Product information

Product name	Test samples	Specifications	Storage conditions
Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)	Human Anterior Nares Swab	1T/kit, 5T/kit, 20T/kit, 40T/kit	4-30 C

Company Profile

Beijing Hotgen Biotech Co., Ltd. (abbreviated as Hotgen Biotech, stock code: 688068) was established in June 2005, which is a high-tech enterprise focusing on the research& development, manufacture and sales of medical and public safety inspection products of in vitro diagnostics (IVD) in the field of biomedicine, as well as landed on the China Sci-Tech innovation board (STAR Market) in September 2019.

After several years of Research& development, Hotgen Biotech has developed an in vitro diagnostic reagent bioactive raw material development platform, a sugar chain abnormal protein detection (sugar capture) R&D technology platform, a Magnetic particles chemiluminescence R&D technology platform, a Up-converting Phosphor R&D technology platform, and a colloidal gold immune layer, The eight major technology platforms, such as the precipitation R&D platform, enzyme-linked immunoassay R&D technology platform, molecular diagnostics R&D platform, and instrument R&D technology platform, form a closed-loop system for in vitro diagnostic R&D and production. Hotgen Biotech has established a complete full level immunodiagnostic technology platform, from high-precision Up-converting Phosphor POCT (UPT series) to small, medium and large single- cartridge chemiluminescence platforms (MQ60 series), and then to large-scale full-automatic chemiluminescence Platform (C2000), which realizes the application of the immune diagnostic platform in the field of full diagnostic scenarios. Supporting products are widely used in the clinical and public safety fields. Specific users include hospitals at all levels, township health centers, third-party testing centers, and medical institutions, as well as medical and health institutions, as well as disease control centers, public security, fire protection, military, ports, food and medicine. Supervision, food and feed enterprises and other public safety fields.

Hotgen Biotech has won the second prize of the National Technology Invention Award, the Gold Medal of Independent Innovation, and the second prize of the Chinese Medical Science and Technology Award; In 2018, Hotgen Biotech was awarded the second prize of the "Technical Invention Category of China Rare Earth Science and Technology Award" by the China Rare Earth Society; Top 100 Private Scientific and Technological Innovations "and" Top 100 Medical Enterprises of the Future "; and" Postdoctoral Scientific Research Workstation "; major science and technology projects in the 12th and 13th five years, 863 plan, science and technology projects of the Beijing Science and Technology Commission, and Zhongguancun High Precision The project's major cutting-edge original technological achievements transformation and industrialization projects.

In the face of the COVID-19 epidemic situation, Beijing Hotgen Biotech Co.,Ltd has organized R&D developed a variety of Covid-19 detection reagents, including Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold), Novel Coronavirus 2019-nCoV Antigen Test (Up-converting Phosphor Immunochromatographic Technology), Coronavirus disease(COVID-19) Antibody Test (Colloidal Gold), Coronavirus disease(COVID-19) IgM/IgG Antibody Rapid Test (Colloidal Gold), Novel Coronavirus 2019-nCoV Antibody Test (Up-converting Phosphor Immunochromatographic Technology), Novel Coronavirus (SARS-CoV-2) Neutralizing Antibodies Test (Colloidal Gold), Coronavirus disease(COVID-19) Nucleic Acid Test Kit (PCR-Fluorescent Probe Method), Disposable virus sampling tube, Nucleic acid Automatic Purification System, Nucleic acid extraction reagent, Biological Sample Releaser kit, etc.It is imperative to fight the epidemic Helping the global fight against epidemics!

Since its establishment, the company has continuously grown its business and has now achieved group development. At present, Hotgen (Langfang), Hotgen (Jilin), Weikekang Technology, Shunjing Biological and many other subsidiaries have been established. Hotgen Biotech marketing and service network has covered all regions of the country. Each province is equipped with professional technical service engineers, academic engineers, etc. who are responsible for pre-sales and after-sales work to meet customer needs. The company takes "developing biotechnology and benefiting human health" as its mission, "quality determines the company's life and death, customers determine the company's success or failure, talents determine the company's rise and fall, innovation determines the company's future" as its core values, and "tests because of me advanced" as its philosophy , High ambitions, technological entrepreneurship, and industrial prosperity!

Declaration of Conformity

Manufacturer:

Name: Beijing Hotgen Biotech Co.,Ltd

Address: 9th building, No.9 Tianfu Street, Biomedical Base,Daxing District, Beijing,
102600, P.R.China

European Representative:

MedNet GmbH

Borkstrasse 10,48163 Muenster,Germany

Product Name:

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Novel Coronavirus 2019-nCoV Antigen Test (Up-converting Phosphor Immunochromatographic Technology)

Classification : **Others of ANNEX II of IVDD**

Conformity Assessment Route: **Annex III**

We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.

General applicable directives:

In Vitro Diagnostic Medical Devices DIRECTIVE 98/79/EC

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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	<p>In Englisch / In English</p> <p>Models A+B: This kit is used for in-vitro qualitative determination of novel coronavirus antigen in human nasal swabs or throat swabs. It is used as rapid investigation for suspected cases of novel coronavirus can also be used as a reconfirmation method for nucleic acid detection in discharged cases.</p> <p>Models C+D (Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) - Saliva): This kit is used for in vitro qualitative determination of novel coronavirus antigen in human saliva. It is used as rapid investigation for suspected cases of novel coronavirus, can also be used as a reconfirmation method for nucleic acid detection in discharged cases.</p> <p>Model E (SARS-CoV-2 Antigen Test (Colloidal Gold): This kit is used for in vitro qualitative determination of SARS-CoV-2 antigen in human anterior nasal swab samples. It is used as rapid investigation for suspected cases of novel coronavirus, can also be used as a reconfirmation method for nucleic acid detection in discharged cases.</p>
--	---

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**Medicines & Healthcare products
Regulatory Agency**

10 South Colonnade
Canary Wharf
London
E14 4PU
United Kingdom

+44 (0) 20 3080 6000
gov.uk/mhra

**Medimap Ltd
2 The Drift
Suffolk
Thurston
IP31 3RT
United Kingdom**

27 February 2021

Dear **Tracey Oakey**

We are pleased to confirm that the application to register or update an existing registration for the following manufacturer, which you submitted on **26 February 2021** has been reviewed:

Application reference: **2021022601194301**

Manufacturer organisation: **Beijing Hotgen Biotech Co Ltd**

Address:

9th Building, No 9 Tianfu Street, Biomedical Base

Daxing District

Beijing

102600

China

Manufacturer registration status: **Registered**

Device(s):

GMDN term	Status	MHRA comment
64787 - SARS-CoV-2 antigen IVD, kit, immunochromatographic test (ICT), rapid	Registered	
64956 - SARS-CoV-2 immunoglobulin A (IgA)/IgG/IgM antibody IVD, kit, immunochromatographic test (ICT), rapid	Registered	
65315 - SARS-CoV-2/Influenza A/B antigen IVD, kit, immunochromatographic test (ICT), rapid	Registered	
65147 - SARS-CoV-2 antigen IVD, kit, chemiluminescent immunoassay	Registered	
65300 - SARS-CoV-2 neutralizing antibody IVD, kit, chemiluminescent immunoassay	Registered	

Please note this letter does not represent any form of accreditation, certification or approval by the UK Competent Authority.

If you stop placing devices on the market or if you are not complying with the Regulations, you should inform us so that we can amend our records. You should be aware that it is an offence to place on the market UKCA or CE marked devices that do not comply with the regulations.

Please inform us of the following chargeable changes:

1. **company/organisation information e.g. name and address**
2. **additional devices (GMDN code or term)**

Please also use the Devices Online Registration Database (DORS) to tell us of the following changes e.g. removal/discontinuation of a device (GMDN) or product from your registration record, change of contact person, telephone number and/or email address, for which payment of our statutory fee does not apply.

Please note that the name and address of manufacturer, UK Responsible Person or Authorised Representative (Northern Ireland only) and devices that have been registered will be published on our [Public Access Registration Database](#) (PARD). This applies to non-in vitro diagnostic devices only.

The account number for your company/organisation is **0000012377**.



Yours sincerely,

Ngozi Onyeukwu
Device registrations service
Devices division
MHRA

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Product Photos

前鼻腔抗原胶体金检测试剂包装信息
Novel Coronavirus 2019-nCoV Antigen Test(Colloidal Gold)
Packing Information

产品名称 Product name	规格/盒 Specifications	单位 Unit	单位包装毛重 Gross weight per unit package
Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)	1T	盒/kit	0.033 kg/盒 0.033 kg / kit

前鼻腔抗原胶体金试剂盒出口包装箱 Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) Export Packing Cartons							
包装箱 / 盒 Packing Carton/ box	长 length cm	宽 Width cm	高 height cm	每箱装盒 数量 Kit quantity per carton	单盒试剂 净重 Net weight of single kit	整箱净重 Net weight of the whole carton	抛重 Throwing weight
纸箱 carton	71	40	39	320盒 320kits	0.033 公斤 0.033 kg	10.56公斤 10.56 kg	18.5-19公 18.5-19 kg

Sensitivity verification of Novel Coronavirus 2019-nCoV

Antigen Test (Colloidal Gold)

Purpose

Use inactivated new coronavirus to evaluate the sensitivity of Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Experimental Materials

1. 1 batch of colloidal gold test paper;
2. Inactivated virus: 10^5 pfu/mL.

Experimental steps

Sample: Mixing ratio of sample diluent

Concentration number	Virus content in sample (pfu/mL)	Sample: Mixing ratio of sample diluent
1	0	1: 9
2	10^2	1: 9
3	2.5×10^2	1: 9
4	5×10^2	1: 9
5	10^3	1: 9
6	10^4	1: 9

1. After mixing the sample and diluent, incubate at room temperature for 1 min.
2. Take 100 μ L of sample and observe the result after 15min reaction.

Test results

Concentration number	Virus content in sample (pfu/mL)	Sample: Mixing ratio of sample diluent	Result
1	0	1: 9	-
2	10^2	1: 9	\pm
3	2.5×10^2	1: 9	+
4	5×10^2	1: 9	+
5	10^3	1: 9	++
6	10^4	1: 9	+++

In conclusion

Colloidal gold experiment results: 10^2 pfu/mL has a shallow band, negative without background, the sensitivity is 2.5×10^2 pfu/mL.

The Key laboratory of Biological Emergency
and Clinical POCT (Beijing)
Aug. 17th, 2020

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Certificate for Safe Transport of Air Cargo

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SINO-Dangerous Goods Transportation Consultant Ltd.

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